

FOR IMMEDIATE RELEASE

Young launches NEW ergonomic hygiene handpiece designed for optimum comfort

Earth City, Mo., Sept. 12, 2014 – Designed to fit the contours of the hygienist’s hand, the new Young Hygiene Handpiece delivers the ergonomic comfort and control clinicians need during polishing.

The first hygiene handpiece ever made by Young Dental is designed with a unique contra-angled connector, which attaches the handpiece to air tubing at a 45-degree angle. This angled connection is designed to help reduce cord drag aggravation and to make the handpiece feel lighter.

Unlike straight-bodied handpieces, the Young Hygiene Handpiece is constructed with a Reverse Radius™ silhouette that’s curved to fit comfortably in the dental hygienist’s hand for a smooth, ergonomic performance. The handpiece’s engraved matte grip rotates 360°, which is designed to help reduce the need to shift the wrist as often during polishing.

“I loved your handpiece for its comfort! The reverse radius nestled in my hand, but most of all, the curved hose attachment and swivel made my wrist feel so much better. What a great design,” said Denise Crowell, RDH, who tested the new handpiece. “The Midwest doesn’t have the curved hose attachment – only the swivel, which causes fatigue. The grip was also superior and I didn’t need to use as much tension.”

The handpiece is designed for optimal performance when used with Young Dental disposable prophylaxis angles, so clinicians can get the most out of their polishing routines. Each handpiece comes with samples of our newest and most compact disposable prophylaxis angle, Vera by Young™.

When it comes to care and maintenance, the Young Hygiene Handpiece nosecone and motor easily detach into two parts for fast preparation and thorough sterilization. The handpiece is proudly made in the USA and is backed by a 2-year warranty with quality guaranteed.

Visit www.youngdental.com for more information about the new Young Hygiene Handpiece.

About Young Dental:

Young Dental is dedicated to providing the dental community with innovative, high-value products that reflect the professionalism of the industry. Based in Earth City, Mo., the company is a leading designer, manufacturer and marketer of consumable supplies, instruments and other products used by dental professionals primarily in preventive dentistry, restorative procedures and instrument sterilization. Young Dental is committed to continuous improvement and the creation of high-quality products.

Visit our website @ www.YoungDental.com

Follow Young Dental on Facebook @ www.facebook.com/YoungDental

Follow Young Dental on Twitter @ www.twitter.com/Young_Dental