



YOUNG™
INFINITY

CORDLESS HYGIENE HANDPIECE



INFINITE CHOICES

INFINITE POSSIBILITIES

Introducing the **NEW** era of **Cordless Hygiene Handpiece**

INFINITE CHOICES

Enjoy the flexibility of a cordless handpiece without sacrificing freedom of choice. Works with any doriot style disposable prophylaxis angle.

EASY ADAPTION

Cordless foot pedal controls speed without restricting range of motion. No need to adjust technique or purchase additional accessories. Unique nosecone design features a 360° swivel, allowing for greater access and maneuverability without compromising form.

EXCEPTIONAL ENGINEERING

Features a light-weight form and well balanced design for ergonomic comfort. Enjoy a full range of motion and greater access without the burdens of cord drag.

INFINITE POSSIBILITIES

Comprehensive kit includes 3 autoclavable nose cones, reducing the need for multiple handpieces per operator. Portable design allows clinicians and practitioners to explore new avenues of patient care outside the traditional boundaries of the operator.

Young™ Infinity Cordless Hygiene System Includes:

1 Cordless Handpiece, 1 Cordless Foot Pedal, 1 Handpiece Cradle, 3 Replacement Nosecones, Charging Cords, 25 Disposable Barrier Sleeves

Item # 295737



YOUNG™ INFINITY

CORDLESS HYGIENE HANDPIECE

“This is a **GAME CHANGER!** Cordless, light weight, easy to use, and I can pick any angle I choose...no restrictions...literally!”

“The Infinity is by far the quietest and smoothest cordless handpiece I have used in practice.”

“This cordless handpiece blew my mind! The **sleek design, power and ergonomic design** were a real plus. And because it can be used with **ANY prophy angle**, I get to use MY favorite angle, not someone else's - now that's #winning!”

“The Young Cordless Handpiece is incredibly **easy to use** and a huge game changer when it comes to cordless handpieces. The speed and power are a hygienist's dream! It is able to power through stain. Being cordless allows for a **quicker, more effective** polish.”



TECHNICAL SPECS:

Weight	118 grams
Dimensions	28mm Dia. x 160mm
Power Source	Rechargeable Lithium Ion Battery
Prophy Angles Accepted	All Doriot Style
Nosecone Swivel	YES
Foot Pedal	YES
Motion	Rotary
Speed Control	Foot Pedal
Speed Range	Variable 500-3000 RPM

OPTIONAL ACCESSORIES:



295721
(1) Handpiece



295791
(1) Removable
Nosecone



295831
(500) Disposable
Barrier Sleeve



295792
(1) Foot Controller
including Stick



295788
(1) Charger Power Cord
and Removable Plugs



295789
(1) Handpiece
Cradle



295790
(3) O-Ring