

INFINITE CHOICES

INFINITE POSSIBILITIES

Introducing the NEW era of **Cordless Hygiene Handpiece**



Works with virtually all disposable prophy angles



Low price point includes 3 autoclavable nosecones and a foot pedal



Quiet operation helps reduce patient anxiety



Simple, effective infection control that complies with CDC guidelines



Battery life lasts nearly **4x LONGER** than the Nupro Freedom Cordless Prophy System¹



YOUNG™
INFINITY

CORDLESS HYGIENE HANDPIECE

“This is a **GAME CHANGER!** Cordless, light weight, easy to use, and I can pick any angle I choose...no restrictions...literally!”

“The Infinity is by far the quietest and smoothest cordless handpiece I have used in practice.”

“This cordless handpiece blew my mind! The **sleek design, power and ergonomic design** were a real plus. And because it can be used with **ANY prophy angle**, I get to use MY favorite angle, not someone else's - now that's #winning!”

“The Young Cordless Handpiece is incredibly **easy to use** and a huge game changer when it comes to cordless handpieces. The speed and power are a hygienist's dream! It is able to power through stain. Being cordless allows for a **quicker, more effective** polish.”



TECHNICAL SPECS:

Weight	118 grams
Dimensions	28mm Dia. x 160mm
Power Source	Rechargeable Lithium Ion Battery
Prophy Angles Accepted	All Doriot Style
Nosecone Swivel	YES
Foot Pedal	YES
Motion	Rotary
Speed Control	Foot Pedal
Speed Range	Variable 500-3000 RPM

OPTIONAL ACCESSORIES:



295721
(1) Handpiece



295791
(1) Removable
Nosecone



295831
(500) Disposable
Barrier Sleeve



295792
(1) Foot Controller
including Stick



295788
(1) Charger Power Cord
and Removable Plugs



295789
(1) Handpiece
Cradle



295790
(3) O-Ring